




## 3. Experiment

Sample:     Cholesterol in CDCl<sub>3</sub> (50mg/ml)  
 Spectrometer: 7 Tesla Magnet with Tecmag HF3 discovery  
 Probe: Nalorac D300-5 OWB 5mm <sup>1</sup>H/<sup>13</sup>C Switchable probe  
<sup>1</sup>H hard pulse:   13.9 kHz (90° = 18 μs @ 5W)  
<sup>1</sup>H decoupling field: 5.6 kHz (90° = 45 μs @ 800 mW)  
<sup>13</sup>C hard pulse:  59.5 kHz (90° = 4.2 μs @ 250 W)  
 τ:   3.5 ms (= 1/(2J<sub>C,H</sub>) = 140 Hz)  
 Spectrum width: ± 6.5kHz  
 Recycling time:   2s  
 Number of scans:  512



## Notes:

1. Before editing the sequence (Fig. 1b), calibrate the 90° pulse widths of <sup>1</sup>H and <sup>13</sup>C using the nutation experiment (see note, "One Pulse Experiment and Pulse Calibration").
2. Set up the WALTZ sequence according to the note, "<sup>13</sup>C NMR Spectra with <sup>1</sup>H WALTZ Decoupling".
3. The center of pulses P2 and P4 (also P3 and P5) should be aligned. Since P2 > P4 (and P3 > P5) P2 (and P3) have to split into 3 pulses. The delay of P2's (and P3's) middle pulse equals to P4 (and P5), and the delay of both sides is (P2 - P4)/2 [and (P3 - P5)/2]. The middle pulse of P2 (and P3) falls on the same event as P4 (and P5).

## 4. Results




Fig. 2. The Editing DEPT<sup>13</sup>C spectra of cholesterol in CDCl<sub>3</sub>. Spectra a, b, and c are correspondent to P3 = 45°, 90°, and 135° in the sequence shown in Fig.1. Spectra d and e are calculated using spectra a, b and c and the equations shown in the figure. The calculation is done by a script, "Editing DEPT".

## 5. References

1. M.R. Bendall, D.M. Doddrell, D.T.Pegg, *J. Am. Chem. Soc.* **1981**, *103*, 4603-4605.
2. D.M. Doddrell, D.T.Pegg, M.R. Bendall, *J. Magn. Res.* **1982**, *48*, 323-327.
3. S. Braun, H.-O. Kalinowski, S. Berger, "150 and More Basic NMR Experiments", Wiley-VCH, 1999, 183-184.