
1. Introduction

The Attached Proton Test (APT) experiment is a common way to assign C-H multiplicities in ¹³C NMR spectra. It provides the information on all sorts of carbons within one experiment. Depending on the number of hydrogens bound to a carbon atom, n, CH_n spin vectors evolve differently after the initial pulse. If the delay τ is set to 1/J, CH and CH₃ vectors have opposite phases compared to C and CH₂. Therefore, the phase of CH and CH₃ peaks is 180° different from C and CH₂ peaks. This example demonstrates the basic procedure of double resonance 1D NMR data acquisition and processing on Tecmag spectrometers.

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
Fig. 1. (a) The ¹³C APT sequence with WALTZ-16 ¹H decoupling. (b) The sequence in the NTNMR sequence editor. The WALTZ-16 sequence is loaded using the asynchronous function.

3. Experiment

Sample:	Cholesterol in CDCl ₃ (50mg/ml)
Spectrometer:	7 Tesla Magnet with Tecmag HF3 discovery
Probe:	Nalocac D300-5 OWB 5mm ¹ H/ ¹³ C Switchable probe
¹ H decoupling field:	5.6 kHz (90° = 45 μs @ 800 mW)
¹³ C hard pulse:	59.5 kHz (90° = 4.2 μs @ 250 W)
τ:	7 ms (= 1/J _{C,H} = 140 Hz)
SW +/-:	± 6.5kHz
Last Delay:	2s
Scans 1D:	512

Notes:

Before editing the sequence (Fig. 1b), you have to calibrate the 90° pulse width of ¹H and ¹³C using the nutation experiment (see notes, "One Pulse Experiment and Pulse Calibration" and set up the WALTZ sequence according to the notes, "¹³C NMR Spectra with ¹H WALTZ Decoupling").

4. Results

Figure 2 shows the ¹³C APT spectrum of cholesterol in CDCl₃ acquired with the sequence introduced in Fig. 1. Note that the signals of CH and CH₃ are positive, but CH₂ and quaternary carbons including the solvent carbon are negative.


Fig. 2. The ¹³C APT spectrum of cholesterol in CDCl₃ using the sequence shown in Fig. 1.

5. References

1. S.L. Patt, D. Ruben and J.N. Shoolery, *J. Magn. Reson.* **1982**, *46*, 535-539.
2. H. J. Jakobsen, O.W. Sorensen, W.S. Brey and P.J. Kanyha, *J. Magn. Reson.* **1982**, *48*, 328.
2. A.M. Torres, T.T. Nakashima, and R.E.D. McClung, *J. Magn. Res. Ser. A* **1993**, *101*, 285-294.
3. S. Braun, H.-O. Kalinowski, S. Berger, "150 and More Basic NMR Experiments", Wiley-VCH, 1999, 165-167.