

1. Introduction


The time constant T_1 is the longitudinal or spin-lattice relaxation time. It is a measure of the time that each nucleus needs to reach equilibrium with its surroundings in the direction of B_0 (z-axis). One popular method to measure T_1 is the inversion-recovery method. The magnetization (M_0) is inverted with a π pulse, then a time τ allowed to relax along the +z-axis. The magnetization, after time τ , is measured applying a $\pi/2$ pulse and then measuring the NMR signal. T_1 is, thus, determined by fitting the magnetization curve or direct calculation of $T_1 = \tau_n / \ln 2$, where τ_n is the τ at node ($M = 0$).

$$M = M_0(1 - 2e^{-\tau/T_1}), \quad (1)$$

2. Pulse sequence

Figure 1 shows an example measuring proton T_1 through $^1\text{H}/^{13}\text{C}$ cross-polarization and ^{13}C observation.

Fig. 1a: the sequence for measuring ^1H T_1 of Gly- $1\text{-}^{13}\text{C}$, ^{15}N . Mixing time is 3ms and rf field is 50 kHz. The decoupling field is less than 50 kHz. b: the actual sequence in NTNMR sequence editor.



(b)

Event Number	1	2	3	4	5	6	7	8	9	10
Name:										
Delay	1u	pd	pi	tau	pw	mix	rd	ad	Acq. Time	Last Delay
F1_Ampl		H90 amp	H90 amp	H90 amp	H90 amp	Hmix amp	Hdec amp	Hdec amp	Hdec amp	
F1_PhMod										
F1_Ph			X	X	X	Y				
F1_Atten		H90 attn	H90 attn	H90 attn	H90 attn	Hmix attn	Hdec attn	Hdec attn	Hdec attn	
F1_TxGate										
F1_PhRst										
Ext_Trig										
F2_Ampl						pw amp				
F2_PhMod										
F2_Ph						π				
F2_Atten						pw attn				
F2_TxGate										
F2_PhRst										
Acq										
Acq_phase										
RX_Blank										
Delay_2D										

0.01s, 0.05s, 0.1s, 0.15s, 0.2s, 0.3s, 0.4s

3. Experiment and results

A set of spectra at different τ 's is obtained in a 2D experiment. A 2D delay table is used in event 4 (see Fig. 1b). Fig. 2 shows the fully labeled Gly-1-¹³C CPMAS spectra acquired with the sequence shown in Fig. 2a and plotted as a function of τ .




Fig. 2 The centerband array of Gly-1-¹³C CPMAS ¹H spectra. Spin speed is 2.5 kHz.

4. T₁ calculation

NTNMR has a built-in data analysis tool which can be used to calculate T₁. First, the spectra are phased such that the first spectrum is inverted. Then the desired peak is highlighted for T₁ calculation. The "NMR Data Analysis" window is open from "Tools" menu, the τ table as "x-axis" and "Intensity" as "source" are selected, the "Draw" button is clicked, and the result will appear in the window as shown in Fig. 3.




Fig. 3. T₁ fitting for glycine-1-¹³C in the "NTNMR Data Analysis" window. Fitting results: T₁ = 0.22 ms.